




# ハードコート・永年保証




### RT-Raman 266nm Laser Longpass Set




- Ex : RET266/3x    ● BS : RT266rdc
- Em: RET270lp    ● 波長帯 : ~ 300nm

### RT-Raman 488nm Laser Longpass Set




- Ex : RET488/4x    ● BS : RT488rdc
- Em: RET493lp    ● 波長帯 : ~ 600nm

### RT-Raman 532nm Laser Longpass Set




- Ex : RET532/4x    ● BS : RT532rdc
- Em: RET537lp    ● 波長帯 : ~ 650nm

### RT-Raman 633nm Laser Longpass Set




- Ex : RET633/5x    ● BS : RT633rdc
- Em: RET638lp    ● 波長帯 : ~ 800nm

### RT-Raman 785nm Laser Longpass Set



- Ex : RET785/6x    ● BS : RT785rdc
- Em: RET792lp    ● 波長帯 : ~ 900nm

### RT-Raman 1064nm Laser Longpass Set



- Ex : RET1064/7x    ● BS : RT1064rdc
- Em: RET1072lp    ● 波長帯 : ~ 1200nm